

AEO Wavelength Meter Specification

Specifications		
Optical parameters	Typical	Note
Wavelength range	350 ~1000 nm	
Accuracy	0.01m	
Resolution	0.005nm	
Sensitivity	20uw	
Power consumption	2 w	
Power supply	10v/60Hz , 220v/50 Hz	
Dimensions	200 (L)x 165 (W) X 70 (H) mm	
Weight	1.2 KG	
Operation temperature	10 °C ~ +3 °C	
Storage temperature:	-10 °C ~ +70 °C	